


mediff®



Differentiating & Evaluation System



Types:


mediff PPC KB-16b	or	mediff SW KB-16b	for blood
mediff PPC KB-30m	or	mediff SW KB-30m	for marrow and blood
mediff PPC KB-30u	or	mediff SW KB-30u	for urine
mediff PPC KB-30x	or	mediff SW KB-30x	for other samples

Specifications:

- ✓ Keyboard with 16 or 30 differentiating plus 6 function keys
- ✓ mediff Software (SW) for Windows CE, XP, Vista, 7
- ✓ LCD colour flat screen 12" with integrated computer (PPC)
- ✓ Differentiating of cells 1x1 to 5x1000 (b, m)
- ✓ Validation out of 1 to 5 counts possible
- ✓ Reference range with graphic display (several possible)
- ✓ Correction option (inserting/deleting of cells & counts)
- ✓ Multiple users possible (for example per count another user)
- ✓ Connection to Laboratory Information System

Options:

- ✓ WBC absolute value distribution with graphic display
- ✓ Morphology evaluation with text (b, m)
- ✓ Display of haematology analyzer results (b, m)
- ✓ Urine data from analyzer (u)
- ✓ Urine analyzer, test stick instrument with uretic parameter, networked (u)
- ✓ Patient data file local or network
- ✓ Patient particulars from LIS
- ✓ Connection to Laboratory Information System (LAN)
- ✓ Enlargement to other probe type
- ✓ Special key lettering
- ✓ Remarks as text
- ✓ Multiple stations (client host)
- ✓ Remote maintenance (Remote Desktop)
- ✓ Barcode system
- ✓ Thermo printer



Advantages:

- ✓ All results on a view
- ✓ Diagrams
- ✓ Optimum validation

30 years innovation for success